

Medical emergency card

I am wearing a required medical device prescribed by my physician.



Medtronic

In an emergency, please contact:

physician name:

phone number:

family member:

phone number:

I have diabetes.

If my behavior is peculiar (appear intoxicated) or I am unconscious, it may be a result of severe low blood sugar.

I am not intoxicated.

Call a medical team for assistance.

If I am able to swallow, give me a source of sugar, (examples: juice, candy, non-diet soft drink).

Healthcare Professionals

I am wearing an external insulin pump.

- The pump is delivering fast-acting insulin at a constant rate.
- To stop insulin delivery, press ACT, use down arrow button to select Suspend, press ACT. Follow instructions on the screen.
- Do not remove the pump battery.
- Do not remove the pump without medical consent.
- If the pump is alarming, follow the instructions on the pump screen or call the local Medtronic Diabetes representative listed below.

Local Medtronic Diabetes Representative

name:

phone number:

Airport information

I am wearing a required medical device prescribed by my physician.



Medtronic

This device is safe for use on commercial airlines. If questioned by airline personnel about use of your Medical Portable Electric Device (M-PED), show them this card.

The transmitter is an M-PED with emission levels that meet RTCA/DO160D Section 21, Category M. Per FAA Advisory Circular #91-21.1B, dated 8/25/06; any medical-portable electronic device (M-PED) that meets RTCA/DO160D Section 21, Category M in any mode “may be used onboard the aircraft without any further testing by the operator”.



This device can withstand exposure to common electrostatic (ESD) and electromagnetic interference (EMI), including airport security systems. During exposure, it can remain connected and powered on.

Do not expose this device to any radiation (such as X-ray, MRI, CT scan).

Patient information:

my name:

my address:

home phone number:

Medical Device Information

medical device type:

device serial number:



Medtronic

Medtronic MiniMed
Northridge, CA 91325
USA
800-646-4633
818-576-5555

E.U. Representative
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
31 (0) 45 566 8000
www.medtronicdiabetes.com